

Anti-Daxx Rabbit pAb

Purified Rabbit Polyclonal Antibody


Catalog # P010452

Product Information

Application	WB, IHC-P/IF (Tissue-P), ICC/IF (Cell), FC, ELISA
Reactivity	Human
Dilution	WB 1:500~1:1,000; IHC-P 1:50~1:100; IF 1:50~1:200; FC 1:50~1:100
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	A synthesized peptide derived from human Daxx
Format	Buffer System: Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Purification: Affinity Chromatography
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-Daxx antibody is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	DAXX, BING2, DAP6, Death domain-associated protein 6, Daxx, hDaxx, ETS1-associated protein 1, EAP1, Fas death domain-associated protein.
Calculated MW	Calculated MW: 81 kDa; Observed MW: 110 kDa
Uniprot ID	Q9UER7
Gene ID	1616
Background	Acts as an adapter protein in a MDM2-DAXX-USP7 complex by regulating the RING-finger E3 ligase MDM2 ubiquitination activity. Under non-stress condition, in association with the deubiquitinating USP7, prevents MDM2 self-ubiquitination and enhances the intrinsic E3 ligase activity of MDM2 towards TP53, thereby promoting TP53 ubiquitination and subsequent proteasomal degradation. Upon DNA damage, its association with MDM2 and USP7 is disrupted.



Western blot analysis of Calreticulin in HeLa lysates using Daxx antibody.