

Anti-ATG7 Rabbit pAb

Purified Rabbit Polyclonal Antibody


Catalog # P013110

Product Information


Application	WB, IHC-P/IF (Tissue-P), ICC/IF (Cell), IP, FC, ELISA
Reactivity	Human
Dilution	WB 1:500~1:1,000; IHC-P 1:50~1:100; IF 1:50~1:200; IP 1:20; FC 1:50~1:100
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	A synthesized peptide derived from human ATG7
Format	Buffer System: Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Purification: Affinity Chromatography
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-ATG7 antibody is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	hAGP7, Ubiquitin-activating enzyme E1-like protein, APG7L.
Calculated MW	Calculated MW: 78 kDa; Observed MW: 78 kDa
Uniprot ID	O95352
Gene ID	10533
Background	The molecular machinery of autophagy was largely discovered in yeast and referred to as autophagy-related (Atg) genes. Formation of the autophagosome involves a ubiquitin-like conjugation system in which Atg12 is covalently bound to Atg5 and targeted to autophagosome vesicles. This conjugation reaction is mediated by the ubiquitin E1-like enzyme Atg7 and the E2-like enzyme Atg10.



Western blot analysis of ATG7 in HepG2 lysates using ATG7 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon using ATG7 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.