

Anti-VAV2 Rabbit pAb

Purified Rabbit Polyclonal Antibody


Catalog # P013022

Product Information


Application	WB, IHC-P/IF (Tissue-P), ICC/IF (Cell), IP, ELISA
Reactivity	Rat, Human, Mouse
Dilution	WB 1:500~1:1,000; IHC-P 1:50~1:100; IF 1:50~1:200; IP 1:20
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	A synthesized peptide derived from human VAV2
Format	Buffer System: Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Purification: Affinity Chromatography
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-VAV2 antibody is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	VAV2, Vav 2 oncogene, VAV-2.
Calculated MW	Calculated MW: 101 kDa; Observed MW: 101 kDa
Uniprot ID	P52735
Gene ID	7410
Background	Vav proteins belong to the Dbl family of guanine nucleotide exchange factors (GEFs) for Rho/Rac small GTPases. The three identified mammalian Vav proteins (Vav1, Vav2 and Vav3) differ in their expression. Vav1 is expressed only in hematopoietic cells and is involved in the formation of the immune synapse. Vav2 and Vav3 are more ubiquitously expressed.



Western blot analysis of VAV2 in 293T lysates using VAV2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using VAV2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.