

Anti-Fas Rabbit pAb

Purified Rabbit Polyclonal Antibody


Catalog # P010852

Product Information

Application	ELISA, WB, IHC-P/IF (Tissue-P), ICC/IF (Cell), FC
Reactivity	Human
Dilution	WB 1:500~1:1,000; IHC-P 1:50~1:100; IF 1:50~1:200; FC 1:50~1:100
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	A synthesized peptide derived from human Fas
Format	Buffer System: Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Purification: Affinity Chromatography
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-Fas antibody is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	FAS, ALPS1A, APO-1, APT1, CD95, FAS1, FASTM, TNFRSF6.
Calculated MW	Calculated MW: 38 kDa; Observed MW: 48 kDa
Uniprot ID	P25445
Gene ID	355
Background	Receptor for TNFSF6/FASLG. The adapter molecule FADD recruits caspase-8 to the activated receptor. The resulting death-inducing signaling complex (DISC) performs caspase-8 proteolytic activation which initiates the subsequent cascade of caspases (aspartate-specific cysteine proteases) mediating apoptosis.



Western blot analysis of Fas in Ramos lysates using Fas antibody.