

Anti-NUR77 Rabbit pAb

Purified Rabbit Polyclonal Antibody


Catalog # P011093

Product Information

Application	WB, ELISA
Reactivity	Human, Mouse, Rat
Dilution	WB 1:500~1:1,000
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	A synthesized peptide derived from human NUR77
Format	Buffer System: Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Purification: Affinity Chromatography
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-NUR77 antibody is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	NR4A1, GFRP1, HMR, NAK1, Nuclear receptor subfamily 4 group A member 1, Early response protein NAK1, Nuclear hormone receptor NUR/77, Nur77, Orphan nuclear receptor HMR, Orphan nuclear receptor TR3, ST-59, Testicular receptor 3.
Calculated MW	Calculated MW: 64 kDa; Observed MW: 64 kDa
Uniprot ID	P22736
Gene ID	3164
Background	Orphan nuclear receptor. May act concomitantly with NURR1 in regulating the expression of delayed-early genes during liver regeneration. Binds the NGFI-B response element (NBRE) 5'-AAAAGGTCA-3' (By similarity). May inhibit NF-kappa-B transactivation of IL2.



Western blot analysis of NUR77 in HepG2 lysates using NUR77 antibody.