

Anti-Wnt2b Rabbit pAb

Purified Rabbit Polyclonal Antibody


Catalog # P013459

Product Information

Application	WB, ICC/IF (Cell), FC, ELISA
Reactivity	Human, Mouse, Rat
Dilution	WB 1:500~1:1,000; IF 1:50~1:200; FC 1:50~1:100
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	A synthesized peptide derived from human Wnt2b
Format	Buffer System: Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Purification: Affinity Chromatography
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-Wnt2b antibody is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	WNT13, Wnt2b, XWNT2.
Calculated MW	Calculated MW: 44 kDa; Observed MW: 44 kDa
Uniprot ID	Q93097
Gene ID	7482
Background	Ligand for members of the frizzled family of seven transmembrane receptors. Probable developmental protein. May be a signaling molecule which affects the development of discrete regions of tissues. Is likely to signal over only few cell diameters.



Western blot analysis of Wnt2b in HepG2 lysates using Wnt2b antibody.