

Anti-NCK1 Rabbit pAb

Purified Rabbit Polyclonal Antibody


Catalog # P012012

Product Information

Application	WB, IHC-P/IF (Tissue-P), ICC/IF (Cell), IP, FC, ELISA
Reactivity	Human, Mouse, Rat
Dilution	WB 1:500~1:1,000; IHC-P 1:50~1:100; IF 1:50~1:200; IP 1:20; FC 1:50~1:100
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	A synthesized peptide derived from human Nck
Format	Buffer System: Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Purification: Affinity Chromatography
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-NCK1 antibody is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	NCK, nck-1, NCKalpha.
Calculated MW	Calculated MW: 43 kDa; Observed MW: 47 kDa
Uniprot ID	P16333
Gene ID	4690
Background	Adapter protein which associates with tyrosine-phosphorylated growth factor receptors or their cellular substrates. Maintains low levels of EIF2S1 phosphorylation by promoting its dephosphorylation by PP1.



Western blot analysis of Nck in HeLa lysates using NCK1 antibody.