

Anti-Bcl6 Rabbit pAb

Purified Rabbit Polyclonal Antibody


Catalog # P010718

Product Information

Application	ICC/IF (Cell), WB, IHC-F/IF (Tissue-F), IHC-P/IF (Tissue-P), ELISA
Reactivity	Human, Mouse, Rat
Dilution	WB 1:500~1:1,000; IHC-P 1:50~1:100; IF 1:50~1:200; ELISA 1:10,000
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human BCL6.
Format	Buffer System: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3. Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-Bcl6 antibody is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	BCL6, BCL5, LAZ3, ZBTB27, ZNF51, B-cell lymphoma 6 protein, BCL-6, B-cell lymphoma 5 protein, BCL-5, Protein LAZ-3, Zinc finger and BTB domain-containing protein 27, Zinc finger protein 51.
Calculated MW	Calculated MW: 79 kDa; Observed MW: 79 kDa
Uniprot ID	P41182
Gene ID	604
Background	Bcl-6, a transcriptional repressor, binds Stat recognition-like DNA elements and influences germinal center development and cell differentiation. Additionally, Bcl-6 negatively regulates NFκB expression, thereby inhibiting NFκB-mediated cellular functions. HDAC- and silent information regulator (SIR)-2-dependent acetylation of Bcl-6 causes downregulation of activity by inhibiting the ability of Bcl-6 to recruit complexes containing histone deacetylases (HDACs).



Western blot analysis of Bcl6 in 293 lysates using Bcl6 antibody.