

Anti-MKNK2 Rabbit pAb

Purified Rabbit Polyclonal Antibody


Catalog # P013690

Product Information

Application	WB, ELISA
Reactivity	Human, Mouse, Rat
Dilution	WB 1:500~1:1,000; ELISA 1:10,000
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	Synthesized peptide derived from the Internal region of human Mnk2.
Format	Buffer System: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3. Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-MKNK2 antibody is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	MKNK2, GPRK7, MNK2, MAP kinase-interacting serine/threonine-protein kinase 2, MAP kinase signal-integrating kinase 2, MAPK signal-integrating kinase 2, Mnk2.
Calculated MW	Calculated MW: 52 kDa; Observed MW: 52 kDa
Uniprot ID	Q9HBH9
Gene ID	2872
Background	Serine/threonine-protein kinase that phosphorylates SFPQ/PSF, HNRNPA1 and EIF4E.



Western blot analysis of MKNK2 in HeLa lysates using MKNK2 antibody.