

Anti-VASP Rabbit pAb

Purified Rabbit Polyclonal Antibody


Catalog # P013465

Product Information


Application	WB, IHC-F/IF (Tissue-F), IHC-P/IF (Tissue-P), ICC/IF (Cell), ELISA
Reactivity	Human, Mouse, Rat, Monkey
Dilution	WB 1:500~1:1,000; IHC-P 1:50~1:100; IF 1:50~1:200; ELISA 1:10,000
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	The antiserum was produced against synthesized peptide derived from human VASP.
Format	Buffer System: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3. Purification: Affinity Chromatography
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-VASP antibody is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	Vasodilator-stimulated phosphoprotein.
Calculated MW	Calculated MW: 40 kDa; Observed MW: 46,50 kDa
Uniprot ID	P50552
Gene ID	7408
Background	VASP (vasodilator-stimulated phosphoprotein), is involved in the maintenance of cytoarchitecture by interacting with actin-like filaments. VASP shares a limited degree of homology with the amino terminus of WASP, which is frequently mutated in WAS patients. An established substrate of cAMP and cGMP dependent kinases, VASP is phosphorylated on a regulatory serine residue 157 and localizes to focal adhesions, microfilaments and highly active regions of the plasma membrane.



Western blot analysis of VASP in various lysates using VASP antibody.



Western blot analysis of VASP in NIH/3T3 lysates treated with forskolin using VASP antibody. The lane on the right is blocked with the synthesized peptide.