

Anti-Phospho-c Myb (Ser11) Rabbit pAb

Purified Rabbit Polyclonal Antibody

Catalog # P012929


Product Information

Application	WB, IHC-P/IF (Tissue-P), ICC/IF (Cell), IP, ELISA
Reactivity	Human
Dilution	WB 1:500~1:1,000; IHC-P 1:50~1:100; IF 1:50~1:200; IP 1:20
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	A synthesized peptide derived from human Myb
Format	Buffer System: Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Purification: Affinity Chromatography
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-Phospho-c Myb (Ser11) antibody is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	MYB, Transcriptional activator Myb, Proto-oncogene c-Myb.
Calculated MW	Calculated MW: 72 kDa; Observed MW: 72,80 kDa
Uniprot ID	P10242
Gene ID	4602
Background	Transcriptional activator; DNA-binding protein that specifically recognize the sequence 5'-YAAC[GT]G-3'. Plays an important role in the control of proliferation and differentiation of hematopoietic progenitor cells.

Validation Images



Western blot analysis of Phospho-vMyb + c-Myb (S11) in Ramos lysates using Phospho-c Myb (Ser11) antibody.