

Anti-Phospho-Raf1 (Ser259) Rabbit pAb

Purified Rabbit Polyclonal Antibody


Catalog # P011150

Product Information

Application	WB, IHC-P/IF (Tissue-P), ELISA
Reactivity	Human, Mouse, Rat
Dilution	WB 1:500~1:1,000; IHC-P 1:50~1:100
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	A synthesized peptide derived from human Phospho-Raf1 (S259)
Format	Buffer System: Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Purification: Affinity Chromatography
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-Phospho-Raf1 (Ser259) antibody is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	RAF1, RAF, RAF proto-oncogene serine/threonine-protein kinase, Proto-oncogene c-RAF, cRaf, Raf-1.
Calculated MW	Calculated MW: 73 kDa; Observed MW: 73 kDa
Uniprot ID	P04049
Gene ID	5894
Background	Raf-1 is a MAP kinase kinase kinase (MAP3K) which functions downstream of the Ras family of membrane associated GTPases to which it binds directly. Once activated Raf-1 can phosphorylate to activate the dual specificity protein kinases MEK1 and MEK2 which in turn phosphorylate to activate the serine/threonine specific protein kinases ERK1 and ERK2.



Western blot analysis of Phospho-Raf1 (S259) in (1) 293T lysates; (2) 293T lysates treated with Okadaic Acid + Calyculin A using Phospho-Raf1 (Ser259) antibody.