

Anti-Phospho-SIRT1 (Ser47) Rabbit pAb

Purified Rabbit Polyclonal Antibody


Catalog # P013112

Product Information

Application	WB, IHC-F/IF (Tissue-F), IHC-P/IF (Tissue-P), ICC/IF (Cell), ELISA
Reactivity	Human
Dilution	WB 1:500~1:1,000; IHC-P 1:50~1:100; IF 1:50~1:200
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	A synthetic phosphopeptide corresponding to residues surrounding Ser47 of human SIRT1
Format	Buffer System: 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-Phospho-SIRT1 (Ser47) antibody is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	SIRT1, SIR2L1, NAD-dependent protein deacetylase sirtuin-1, hSIRT1, Regulatory protein SIR2 homolog 1, SIR2-like protein 1, hSIR2.
Calculated MW	Calculated MW: 82 kDa; Observed MW: 120 kDa
Uniprot ID	Q96EB6
Gene ID	23411
Background	This gene encodes a member of the sirtuin family of proteins, homologs to the yeast Sir2 protein. Members of the sirtuin family are characterized by a sirtuin core domain and grouped into four classes. The functions of human sirtuins have not yet been determined; however, yeast sirtuin proteins are known to regulate epigenetic gene silencing and suppress recombination of rDNA.



Western blot analysis of Phospho-SIRT1 (Ser47) in A549 lysates using Phospho-SIRT1 (Ser47) antibody.