

Anti-p38 Mouse mAb

Purified Mouse Monoclonal Antibody


Catalog # M013054

Product Information

Application	IHC-P/IF (Tissue-P), ELISA
Reactivity	Human, Rat, Mouse (Cell)
Dilution	IHC-P 1:50~1:100
Host	Mouse
Clonality	Monoclonal
Clone No.	50L90L82
Isotype	IgG1
Label	Unconjugated
Immunogen	Synthetic peptide conjugated to KLH.
Format	Buffer System: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3. Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-p38 antibody [50L90L82] is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	MAPK14, CSBP, CSBP1, CSBP2, CSPB1, MXI2, SAPK2A, Mitogen-activated protein kinase 14, MAP kinase 14, MAPK 14, Cytokine suppressive anti-inflammatory drug-binding protein, CSAID-binding protein, CSBP, MAP kinase MXI2, MAX-interacting protein.
Uniprot ID	Q16539
Gene ID	1432
Background	Responds to activation by environmental stress, pro-inflammatory cytokines and lipopolysaccharide (LPS) by phosphorylating a number of transcription factors, such as ELK1 and ATF2 and several downstream kinases, such as MAPKAPK2 and MAPKAPK5. Plays a critical role in the production of some cytokines, for example IL-6.



Immunohistochemistry analysis of paraffin-embedded Human Lung Carcinoma Tissue using p38 antibody [50L90L82]. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.