

Anti-NFAT5 Rabbit pAb

Purified Rabbit Polyclonal Antibody


Catalog # P013350

Product Information


Application	WB, IHC-P/IF (Tissue-P), ELISA
Reactivity	Human, Rat, Mouse
Dilution	WB 1:500~1:1,000; IHC-P 1:50~1:100; ELISA 1:10,000
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	Synthetic peptide from human protein.
Format	Buffer System: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3. Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-NFAT5 antibody is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information


Synonyms	KIAA0827 TONEBP.
Calculated MW	Calculated MW: 166 kDa; Observed MW: 150-160 kDa
Uniprot ID	O94916
Gene ID	10725
Background	Transcription factor involved, among others, in the transcriptional regulation of osmoprotective and inflammatory genes.



Western blot analysis of NFAT5 in 823 , 293TUV, HELA lysates using NFAT5 antibody.



Immunohistochemistry analysis of paraffin-embedded Human brain using NFAT5 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human brain using NFAT5 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.