

## Anti-C1RL Rabbit pAb

Purified Rabbit Polyclonal Antibody


Catalog # P013359

### Product Information

|             |                                                                                                                                 |
|-------------|---------------------------------------------------------------------------------------------------------------------------------|
| Application | WB, ELISA                                                                                                                       |
| Reactivity  | Human                                                                                                                           |
| Dilution    | WB 1:500~1:1,000; ELISA 1:10,000                                                                                                |
| Host        | Rabbit                                                                                                                          |
| Clonality   | Polyclonal                                                                                                                      |
| Isotype     | IgG                                                                                                                             |
| Label       | Unconjugated                                                                                                                    |
| Immunogen   | Synthesized peptide derived from human C1RL.                                                                                    |
| Format      | Buffer System: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3. Purification: Affinity Purified. |
| Storage     | Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.     |
| Precautions | Anti-C1RL antibody is for research use only and not for use in diagnostic or therapeutic procedures.                            |

### Protein Information

|               |                                                                                                                                   |
|---------------|-----------------------------------------------------------------------------------------------------------------------------------|
| Synonyms      | Complement C1r subcomponent-like protein (C1r-LP, C1r-like protein, EC 3.4.21.-, C1r-like serine protease analog protein, CLSPa). |
| Calculated MW | Calculated MW: 53 kDa; Observed MW: 53 kDa                                                                                        |
| Uniprot ID    | Q9NZP8                                                                                                                            |
| Gene ID       | 51279                                                                                                                             |
| Background    | Mediates the proteolytic cleavage of HP/haptoglobin in the endoplasmic reticulum.                                                 |



Western blot analysis of C1RL in hela lysates using C1RL antibody.