

Anti-Cleaved-Caspase 9 Rabbit pAb

Purified Rabbit Polyclonal Antibody


Catalog # P012889

Product Information

Application	ELISA, ICC/IF (Cell), WB, IHC-F/IF (Tissue-F), IHC-P/IF (Tissue-P)
Reactivity	Mouse, Rat
Dilution	WB 1:500~1:1,000; IHC-P 1:50~1:100; IF 1:50~1:200; ELISA 1:10,000
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	The antiserum was produced against synthesized peptide derived from rat Caspase 9.
Format	Buffer System: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3. Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-Cleaved-Caspase 9 antibody is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	CASP9, MCH6, Caspase-9, CASP-9, Apoptotic protease Mch-6, Apoptotic protease-activating factor 3, APAF-3, ICE-like apoptotic protease 6, ICE-LAP6.
Calculated MW	Calculated MW: 50,17 kDa; Observed MW: 17 kDa
Uniprot ID	Q9JHK1
Gene ID	58918
Background	Involved in the activation cascade of caspases responsible for apoptosis execution. Binding of caspase-9 to Apaf-1 leads to activation of the protease which then cleaves and activates caspase-3.



Western blot analysis of Cleaved-Caspase 9 in various lysates using Cleaved-Caspase 9 antibody.