

Anti-Phospho-beta Catenin (Ser37) Rabbit pAb

Purified Rabbit Polyclonal Antibody


Catalog # P013101

Product Information

Application	WB, IHC-F/IF (Tissue-F), IHC-P/IF (Tissue-P), ICC/IF (Cell), ELISA
Reactivity	Human, Mouse, Rat
Dilution	WB 1:500~1:1,000; IHC-P 1:50~1:100; IF 1:50~1:200; ELISA 1:10,000
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	The antiserum was produced against synthesized peptide derived from human Catenin-beta around the phosphorylation site of Ser37.
Format	Buffer System: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3. Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-Phospho-beta Catenin (Ser37) antibody is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	CTNNB1, CTNNB, OK/SW-cl35, Catenin beta-1, Beta-catenin.
Calculated MW	Calculated MW: 85 kDa; Observed MW: 85 kDa
Uniprot ID	P35222
Gene ID	1499
Background	Beta-catenin is an adherens junction protein. Adherens junctions (AJs; also called the zonula adherens) are critical for the establishment and maintenance of epithelial layers, such as those lining organ surfaces. AJs mediate adhesion between cells, communicate a signal that neighboring cells are present, and anchor the actin cytoskeleton. In serving these roles, AJs regulate normal cell growth and behavior.



Western blot analysis of Phospho-beta Catenin (Ser37) in various lysates using Phospho-beta Catenin (Ser37) antibody.