

Anti-Phospho-DDIT3 (Ser30) Rabbit pAb

Purified Rabbit Polyclonal Antibody


Catalog # P010277

Product Information

Application	WB, ELISA
Reactivity	Human
Dilution	WB 1:500~1:1,000
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	Synthetic peptide of human DDIT3
Format	Buffer System: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3. Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-Phospho-DDIT3 (Ser30) antibody is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	CHOP, CHOP10, GADD153, DNA damage-inducible transcript 3 protein, DDIT-3, C/EBP-homologous protein, CHOP, C/EBP-homologous protein 10, CHOP-10, Growth arrest and DNA damage-inducible protein GADD153.
Calculated MW	Calculated MW: 19 kDa; Observed MW: 19 kDa
Uniprot ID	P35638
Gene ID	1649
Background	Inhibits the DNA-binding activity of C/EBP and LAP by forming heterodimers that cannot bind DNA.



Western blot analysis of Phospho-DDIT3 (Ser30) in Jurkat lysates using Phospho-DDIT3 (Ser30) antibody.