

Anti-PKC delta Rabbit pAb

Purified Rabbit Polyclonal Antibody


Catalog # P012290

Product Information

Application	WB, IHC-P/IF (Tissue-P), ICC/IF (Cell), ELISA
Reactivity	Human, Mouse, Rat
Dilution	WB 1:500~1:1,000; IHC-P 1:50~1:100; IF 1:50~1:200
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	A synthesized peptide derived from human PKC delta
Format	Buffer System: Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Purification: Affinity Chromatography
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-PKC delta antibody is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	Protein kinase C delta type, Tyrosine-protein kinase PRKCD, nPKC-delta.
Calculated MW	Calculated MW: 78 kDa; Observed MW: 78 kDa
Uniprot ID	Q05655
Gene ID	5580
Background	Protein kinase C (PKC) is a family of serine- and threonine-specific protein kinases that can be activated by calcium and the second messenger diacylglycerol. PKC family members phosphorylate a wide variety of protein targets and are known to be involved in diverse cellular signaling pathways.



Western blot analysis of PKC delta in HeLa lysates using PKC delta antibody.