

Anti-GRB2 Rabbit pAb

Purified Rabbit Polyclonal Antibody


Catalog # P011587

Product Information


Application	WB, IHC-P/IF (Tissue-P), IP, ELISA
Reactivity	Human, Mouse, Rat
Dilution	WB 1:500~1:1,000; IHC-P 1:50~1:100; IP 1:20
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	A synthetic peptide of human GRB2
Format	Buffer System: 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-GRB2 antibody is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information


Synonyms	GRB2, ASH, Growth factor receptor-bound protein 2, Adapter protein GRB2, Protein Ash, SH2/SH3 adapter GRB2.
Calculated MW	Calculated MW: 25 kDa; Observed MW: 25 kDa
Uniprot ID	P62993
Gene ID	2885
Background	Adapter protein that provides a critical link between cell surface growth factor receptors and the Ras signaling pathway.



Western blot analysis of GRB2 in PC-12, C6 lysates using GRB2 antibody.



Western blot analysis of GRB2 in MCF-7, 3T3, Raw264.7, HeLa lysates using GRB2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using GRB2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.