

Anti-JNK1 Rabbit mAb

Purified Rabbit Monoclonal Antibody


Catalog # R011131

Product Information

Application	WB, ELISA
Reactivity	Human, Rat
Dilution	WB 1:500~1:1,000
Host	Rabbit
Clonality	Monoclonal
Clone No.	54M05K37
Isotype	IgG
Label	Unconjugated
Immunogen	Recombinant protein of mouse JNK1
Format	Buffer System: 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-JNK1 antibody [54M05K37] is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	JNK, JNK1, PRKM8, SAPK1, JNK-46, JNK1A2, SAPK1c, JNK21B1/2.
Calculated MW	Calculated MW: 44 kDa; Observed MW: 46,54 kDa
Uniprot ID	Q91Y86
Gene ID	26419
Background	Serine/threonine-protein kinase involved in various processes such as cell proliferation, differentiation, migration, transformation and programmed cell death.



Western blot analysis of JNK1 in Jurkat, rat Brain, C6 lysates using JNK1 antibody.