

Anti-Phospho-FAK (Tyr397) Rabbit mAb

Purified Rabbit Monoclonal Antibody


Catalog # R012394

Product Information

Application	WB, ELISA
Reactivity	Rat
Dilution	WB 1:500~1:1,000
Host	Rabbit
Clonality	Monoclonal
Clone No.	23L06L21
Isotype	IgG
Label	Unconjugated
Immunogen	A synthetic phosphopeptide corresponding to residues surrounding Tyr397 of human FAK
Format	Buffer System: 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-Phospho-FAK (Tyr397) antibody [23L06L21] is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	PTK2, FAK, FAK1, Focal adhesion kinase 1, FADK 1, Focal adhesion kinase-related nonkinase, FRNK, Protein phosphatase 1 regulatory subunit 71, PPP1R71, Protein-tyrosine kinase 2, p125FAK, pp125FAK.
Calculated MW	Calculated MW: 119 kDa; Observed MW: 119 kDa
Uniprot ID	Q05397
Gene ID	5747
Background	This gene encodes a cytoplasmic protein tyrosine kinase which is found concentrated in the focal adhesions that form between cells growing in the presence of extracellular matrix constituents. The encoded protein is a member of the FAK subfamily of protein tyrosine kinases but lacks significant sequence similarity to kinases from other subfamilies. Activation of this gene may be an important early step in cell growth and intracellular signal transduction pathways triggered in response to certain neural peptides or to cell interactions with the extracellular matrix.



Western blot analysis of Phospho-FAK (Tyr397) in C6 lysates using Phospho-FAK (Tyr397) antibody.