

Anti-IKB alpha Rabbit mAb

Purified Rabbit Monoclonal Antibody


Catalog # R011223

Product Information


Application	IP, ELISA, WB, IHC-P/IF (Tissue-P)
Reactivity	Mouse, Human, Rat
Dilution	WB 1:500~1:1,000; IHC-P 1:50~1:100; IP 1:20
Host	Rabbit
Clonality	Monoclonal
Clone No.	46K76M90
Isotype	IgG
Label	Unconjugated
Immunogen	A synthetic peptide of human IKB alpha
Format	Buffer System: 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-IKB alpha antibody [46K76M90] is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information


Synonyms	NFKBIA, IKBA, MAD3, NFKBI, NF-kappa-B inhibitor alpha, I-kappa-B-alpha, IkB-alpha, IkappaBalphalpha, Major histocompatibility complex enhancer-binding protein MAD3.
Calculated MW	Calculated MW: 36 kDa; Observed MW: 39 kDa
Uniprot ID	P25963
Gene ID	4792
Background	NFKBI (MIM 164011) or NFKB2 (MIM 164012) is bound to REL (MIM 164910), RELA (MIM 164014), or RELB (MIM 604758) to form the NFKB complex. The NFKB complex is inhibited by I-kappa-B proteins (NFKBIA or NFKBIB, MIM 604495), which inactivate NF-kappa-B by trapping it in the cytoplasm.




Western blot analysis of IKB alpha in PC-12, C6 lysates using IKB alpha antibody.



Western blot analysis of IKB alpha in Raw264.7 lysates using IKB alpha antibody.



Western blot analysis of IKB alpha in 3T3, HeLa lysates using IKB alpha antibody



Immunohistochemistry analysis of paraffin-embedded Human tonsil using IKB alpha antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.