

Anti-Cyclin E2 Rabbit mAb

Purified Recombinant Rabbit Monoclonal Antibody


Catalog # R010460

Product Information


Application	WB, IHC-P/IF (Tissue-P), IF (Cell)/ICC, ELISA
Reactivity	Human
Dilution	WB 1:1,000~1:3,000; IHC-P 1:100~1:200; IF 1:100
Host	Rabbit
Clonality	Monoclonal
Clone No.	52K21M05
Isotype	IgG
Label	Unconjugated
Immunogen	A synthesized peptide derived from human Cyclin E2
Format	Affinity purified monoclonal antibody supplied in PBS with 0.01% sodium azide and 50% glycerol, pH 7.3.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-Cyclin E2 Rabbit mAb [52K21M05] is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information


Synonyms	CCN E2, CCNE2, CCNE2 protein, CYC E2, CYCE 2, CYCE2, CyclinE2, G1/S specific cyclin E2.
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Uniprot ID	O96020
Gene ID	9134
Background	The protein encoded by this gene belongs to the highly conserved cyclin family, whose members are characterized by a dramatic periodicity in protein abundance through the cell cycle. Cyclins function as regulators of CDK kinases. Different cyclins exhibit distinct expression and degradation patterns which contribute to the temporal coordination of each mitotic event. This cyclin forms a complex with and functions as a regulatory subunit of CDK2. This cyclin has been shown to specifically interact with CIP/KIP family of CDK inhibitors, and plays a role in cell cycle G1/S transition. The expression of this gene peaks at the G1-S phase and exhibits a pattern of tissue specificity distinct from that of cyclin E1. A significantly increased expression level of this gene was observed in tumor-derived cells. [provided by RefSeq, Jul 2008]



Western Blot - Anti-Cyclin E2 Rabbit mAb [52K21M05]
All lanes: R010460 at 1:3,000 dilution
Lane 1: HeLa (Human cervix adenocarcinoma epithelial cell) whole cell lysates
Lane 2: HCT116 (Human colorectal carcinoma epithelial cell) whole cell lysates
Lane 3: A431 (Human epidermoid teratoma cell line) whole cell lysates
Lane 4: T24 (Human bladder cancer epithelial cell) whole cell lysates
Lysates/proteins at 10 µg per lane.
Secondary antibody: Goat Anti-Rabbit IgG(H+L), HRP Conjugated (Cat. No. LF102) at 1:5,000 dilution
Predicted band size: 47 kDa
Observed band size: 47 kDa
Developed using the ECL technique (Cat. No. SQ201).



Immunohistochemistry - Anti-Cyclin E2 Rabbit mAb [52K21M05]
Sample: Paraformaldehyde-fixed, paraffin embedded human breast carcinoma tissue
Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0) for 30 mins.
Primary antibody: R010460 at 1:200 dilution
Secondary antibody: Goat Anti-Rabbit IgG (H+L), HRP conjugated at 1:1,000 dilution
DAB was used as the chromogen.
Counter stained with hematoxylin.
Positive/negative staining were presented.
Only the secondary antibody was used as the negative control.



Immunofluorescence - Anti-Cyclin E2 Rabbit mAb [52K21M05]
Sample: HeLa cells
The cells were fixed with 4% paraformaldehyde (10 min), permeabilized with 0.5% Triton X-100 for 10 minutes and then blocked with 5% BSA in 0.1% PBS-Tween for 0.5 hours.
Primary antibodies: R010460 at 1:100 dilution and α-tubulin Mouse Monoclonal Antibody (Cat. No. LF209) at 1:100 dilution
Secondary antibodies: Goat anti-Rabbit (488) at 1:1,000 dilution (shown in green) and Goat anti-Mouse (555) at 1:1,000 dilution (shown in red)
Nuclei were stained with DAPI (shown in blue).