

Anti-FADD Rabbit mAb

Purified Rabbit Monoclonal Antibody


Catalog # R012706

Product Information

Application	WB, IP, ELISA
Reactivity	Mouse, Rat
Dilution	WB 1:500~1:1,000; IP 1:20
Host	Rabbit
Clonality	Monoclonal
Clone No.	89L64L96
Isotype	IgG
Label	Unconjugated
Immunogen	A synthetic peptide of human FADD
Format	Buffer System: 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-FADD antibody [89L64L96] is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	FADD, MORT1, GIG3, Protein FADD, FAS-associated death domain protein, FAS-associating death domain-containing protein, Growth-inhibiting gene 3 protein, Mediator of receptor induced toxicity.
Calculated MW	Calculated MW: 23 kDa; Observed MW: 28 kDa
Uniprot ID	Q13158
Gene ID	8772
Background	Apoptotic adaptor molecule that recruits caspase-8 or caspase-10 to the activated Fas (CD95) or TNFR-1 receptors. The resulting aggregate called the death-inducing signaling complex (DISC) performs caspase-8 proteolytic activation. Active caspase-8 initiates the subsequent cascade of caspases mediating apoptosis.



Western blot analysis of FADD in rat Brain, 3T3 lysates using FADD antibody.