

Anti-Phospho-LAT (Tyr200) Rabbit mAb

Purified Rabbit Monoclonal Antibody


Catalog # R010216

Product Information

Application	WB, ELISA
Reactivity	Human
Dilution	WB 1:500~1:1,000
Host	Rabbit
Clonality	Monoclonal
Clone No.	96L72M05
Isotype	IgG
Label	Unconjugated
Immunogen	A synthetic phosphopeptide corresponding to residues surrounding Tyr200 of human LAT
Format	Buffer System: 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-Phospho-LAT (Tyr200) antibody [96L72M05] is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	LAT, Linker for activation of T-cells family member 1, 36 kDa phospho-tyrosine adapter protein, pp36, p36-38.
Calculated MW	Calculated MW: 28 kDa; Observed MW: 36-38 kDa
Uniprot ID	O43561
Gene ID	27040
Background	Required for TCR (T-cell antigen receptor)- and pre-TCR-mediated signaling, both in mature T-cells and during their development. Involved in FCGR3 (low affinity immunoglobulin gamma Fc region receptor III)-mediated signaling in natural killer cells and FCER1 (high affinity immunoglobulin epsilon receptor)-mediated signaling in mast cells. Couples activation of these receptors and their associated kinases with distal intracellular events such as mobilization of intracellular calcium stores, PKC activation, MAPK activation or cytoskeletal reorganization through the recruitment of PLCG1, GRB2, GRAP2, and other signaling molecules.



Western blot analysis of Phospho-LAT (Tyr200) in Jurkat lysates using Phospho-LAT (Tyr200) antibody [96L72M05].