

Anti-NUMB Rabbit mAb

Purified Recombinant Rabbit Monoclonal Antibody


Catalog # R012880

Product Information


Application	WB, IHC-P/IF (Tissue-P), IF (Cell)/ICC, ELISA
Reactivity	Rat, Human, Mouse
Dilution	WB 1:1,000~1:2,000; IHC-P 1:100~1:200; IF 1:100~1:200
Host	Rabbit
Clonality	Monoclonal
Clone No.	64L17L98
Isotype	IgG
Label	Unconjugated
Immunogen	Recombinant protein of human NUMB
Format	Affinity purified monoclonal antibody supplied in PBS with 0.01% sodium azide and 50% glycerol, pH 7.3.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-NUMB Rabbit mAb [64L17L98] is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information


Synonyms	c14_5527, C14orf41, FLJ31314, h Numb, h-Numb, Numb, Numb homolog (Drosophila), Numb homolog, Numb protein homolog, NUMB_HUMAN, Protein numb homolog, Protein S171, S171.
Calculated MW	Calculated MW: 66, 72 kDa; Observed MW: 66, 72 kDa
Uniprot ID	P49757
Gene ID	8650
Background	Numb acts as a negative regulator of Notch signaling by promoting ubiquitination and degradation of Notch. The protein is asymmetrically segregated into one daughter cell during cell division, producing two daughter cells with different responses to Notch signaling and different cell fates. The localization of Numb can also be regulated by G-protein coupled receptor (GPCR) and PKC signaling.
Cellular Location	Membrane.




Western Blot - Anti-NUMB Rabbit mAb [64L17L98]
All lanes: R012880 at 1:1,000 dilution
Lane 1: HeLa (human cervix adenocarcinoma epithelial cell) whole cell lysates
Lane 2: Jurkat (human T lymphocytic leukemia cell) whole cell lysates
Lane 3: HCT116 (human colorectal carcinoma epithelial cell) whole cell lysates
Lane 4: SCC-9 (human tongue squamous carcinoma epithelial cell) whole cell lysates
Lane 5: U2OS (human osteosarcoma epithelial cell) whole cell lysates
Lane 6: SW620 (human colorectal carcinoma epithelial cell) whole cell lysates
Lane 7: 293T (human embryonic kidney cell) whole cell lysates
Lane 8: HepG2 (human hepatocellular carcinoma epithelial cell) whole cell lysates
Lysates/proteins at 10 µg per lane.
Secondary antibody: Goat Anti-Rabbit IgG(H+L), HRP Conjugated (Cat. No. LF102) at 1:5,000 dilution
Predicted band size: 66 kDa, 72 kDa
Observed band size: 66 kDa, 72 kDa
Developed using the ECL technique (Cat. No. SQ201).




Immunofluorescence - Anti-NUMB Rabbit mAb [64L17L98]
Sample: HeLa cells
The cells were fixed with 4% paraformaldehyde (10 min), permeabilized with 0.5% Triton X-100 for 10 minutes and then blocked with 5% BSA in 0.1% PBS-Tween for 0.5 hours.
Primary antibodies: R012880 at 1:100 dilution and α -tubulin Mouse Monoclonal Antibody (Cat. No. LF209) at 1:100 dilution
Secondary antibodies: Goat anti-Rabbit (488) at 1:1,000 dilution (shown in green) and Goat anti-Mouse (555) at 1:1,000 dilution (shown in red)
Nuclei were stained with DAPI (shown in blue).



Immunofluorescence - Anti-NUMB Rabbit mAb [64L17L98]
Sample: A431 cells
The cells were fixed with 4% paraformaldehyde (10 min), permeabilized with 0.5% Triton X-100 for 10 minutes and then blocked with 5% BSA in 0.1% PBS-Tween for 0.5 hours.
Primary antibodies: R012880 at 1:100 dilution
Secondary antibodies: Goat anti-Rabbit (488) at 1:1,000 dilution (shown in green)
Nuclei were stained with DAPI (shown in blue).



Immunohistochemistry - Anti-NUMB Rabbit mAb [64L17L98]
Sample: Paraformaldehyde-fixed, paraffin embedded human colorectal carcinoma tissue
Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0) for 30 mins.
Primary antibody: R012880 at 1:200 dilution
Secondary antibody: Goat Anti-Rabbit IgG (H+L), HRP conjugated at 1:1,000 dilution
DAB was used as the chromogen.
Counter stained with hematoxylin.
Positive/negative staining were presented.
Only the secondary antibody was used as the negative control.



Immunohistochemistry - Anti-NUMB Rabbit mAb [64L17L98]

Sample: Paraformaldehyde-fixed, paraffin embedded mouse kidney tissue

Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0) for 30 mins.

Primary antibody: R012880 at 1:200 dilution


Secondary antibody: Goat Anti-Rabbit IgG (H+L), HRP conjugated at 1:1,000 dilution

DAB was used as the chromogen.

Counter stained with hematoxylin.

Positive/negative staining were presented.

Only the secondary antibody was used as the negative control.



Western Blot - Anti-NUMB Rabbit mAb [64L17L98]

All lanes: R012880 at 1:1,000 dilution

Lane 1: C2C12 (Mouse myoblasts epithelial cell) whole cell lysates

Lane 2: Raw264.7 (Mouse mononuclear macrophage leukemia cell) whole cell lysates

Lane 3: PC-12 (Rat adrenal pheochromocytoma epithelial cell) whole cell lysates

Lysates/proteins at 10 µg per lane.

Secondary antibody: Goat Anti-Rabbit IgG(H+L), HRP Conjugated (Cat. No. LF102) at 1:5,000 dilution

Predicted band size: 66, 72 kDa

Observed band size: 66, 72 kDa

Developed using the ECL technique (Cat. No. SQ201).