

Anti-Olfactory Marker Protein Rabbit mAb

Purified Rabbit Monoclonal Antibody


Catalog # R012739

Product Information

Application	WB, IP, ELISA
Reactivity	Mouse, Rat
Dilution	WB 1:500~1:1,000; IP 1:20
Host	Rabbit
Clonality	Monoclonal
Clone No.	23L70K02
Isotype	IgG
Label	Unconjugated
Immunogen	Recombinant protein of mouse Olfactory Marker Protein
Format	Buffer System: 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA Purification: Affinity Purified.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-Olfactory Marker Protein antibody [23L70K02] is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	Olfactory marker protein.
Calculated MW	Calculated MW: 19 kDa; Observed MW: 19 kDa
Uniprot ID	Q64288
Gene ID	18378
Background	May act as a modulator of the olfactory signal-transduction cascade.



Western blot analysis of Olfactory Marker Protein in rat Brain, MCF-7 lysates using Olfactory Marker Protein antibody.