

Anti-FLI1 Rabbit mAb

Purified Recombinant Rabbit Monoclonal Antibody


Catalog # R014174

Product Information

Application	WB, ELISA
Reactivity	Human, Mouse
Dilution	WB 1:1,000~1:2,000
Host	Rabbit
Clonality	Monoclonal
Clone No.	93I67I14
Isotype	IgG
Label	Unconjugated
Immunogen	A synthesized peptide derived from human FLI1
Format	Affinity purified monoclonal antibody supplied in PBS with 0.01% sodium azide and 50% glycerol, pH 7.3.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-FLI1 Rabbit mAb [93I67I14] is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	ERGB transcription factor, Ewing Sarcoma breakpoint region 2, EWSR 2, EWSR2, FLI 1, FLI 1 proto oncogene, Fli-1 proto-oncogene, ETS transcription factor, FLI1, FLI1 EWS fusion gene, FLI1 proto oncogene, FLI1_HUMAN, Friend leukemia integration 1 (FLI1) transcription factor, Friend leukemia integration 1 transcription factor, Friend leukemia virus integration 1, Proto-oncogene Fli-1, Sarcoma breakpoint region 2 (EWSR2), SIC 1, SIC1, Transcription factor ERGB, Viral integration region FLI1, Viral integration region FLI1, mouse, homolog of.
Calculated MW	Calculated MW: 51 kDa; Observed MW: 51 kDa
Uniprot ID	Q01543
Gene ID	2313
Background	FLI1 Sequence-specific transcriptional activator. Recognizes the DNA sequence 5'-C[CA]GGAAGT-3'. Belongs to the ETS family. Can form homodimers or heterodimers with ETV6/TEL1. 2 isoforms of the human protein are produced by alternative splicing.
Cellular Location	Nucleus.



Western Blot - Anti-FLI1 Rabbit mAb [93167114]

All lanes: R014174 at 1:1,000 dilution

Lane 1: HeLa (Human cervix adenocarcinoma epithelial cell) whole cell lysates

Lane 2: HepG2 (Human hepatocarcinoma epithelial cell) whole cell lysates

Lane 3: C2C12 (Mouse myoblasts epithelial cell) whole cell lysates

Lysates/proteins at 10 µg per lane.

Secondary antibody: Goat Anti-Rabbit IgG(H+L), HRP Conjugated (Cat. No. LF102) at 1:5,000 dilution

Predicted band size: 51 kDa

Observed band size: 51 kDa

Developed using the ECL technique (Cat. No. SQ201).