

Anti-Phospho-SQSTM1/p62 (Ser349) Rabbit pAb

Purified Rabbit Polyclonal Antibody


Catalog # P104621

Product Information

Application	WB, ELISA
Reactivity	Mouse
Dilution	WB 1:100~1:500
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Label	Unconjugated
Immunogen	A synthetic phosphorylated peptide around S349 of human SQSTM1/p62 (NP_003891.1).
Format	Affinity purified polyclonal antibody supplied in PBS with 0.02% sodium azide and 50% glycerol, pH 7.3.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-Phospho-SQSTM1/p62 (Ser349) Rabbit pAb is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information

Synonyms	p60; p62; A170; DMRV; OSIL; PDB3; ZIP3; p62B; NADGP; FTDALS3; Phospho-SQSTM1/p62-S349.
Calculated MW	Calculated MW: 48 kDa; Observed MW: 62 kDa
Uniprot ID	Q13501
Gene ID	8878
Background	This gene encodes a multifunctional protein that binds ubiquitin and regulates activation of the nuclear factor kappa-B (NF-κB) signaling pathway. The protein functions as a scaffolding/adaptor protein in concert with TNF receptor-associated factor 6 to mediate activation of NF-κB in response to upstream signals. Alternatively spliced transcript variants encoding either the same or different isoforms have been identified for this gene. Mutations in this gene result in sporadic and familial Paget disease of bone.



Western blot analysis of lysates from NIH/3T3 cells, using Phospho-SQSTM1/p62-S349 Rabbit pAb (P104621) at 1:500 dilution. NIH/3T3 cells were treated by MG132(50 μ M) at 37°C for 90 minutes.


Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (LF102) at 1:10,000 dilution.

Lysates/proteins: 25 μ g per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Kit.

Exposure time: 90s.



Western blot analysis of various lysates using Phospho-SQSTM1/p62-S349 Rabbit pAb (P104621) at 1:500 dilution. HeLa cells and C6 cells were treated by MG132(50 μ M) at 37°C for 90 minutes.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (LF102) at 1:10,000 dilution.

Lysates/proteins: 25 μ g per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Enhanced Kit.

Exposure time: 90s.