

PE Mouse IgG2b, κ Isotype Ctrl Antibody

Purified PE-conjugated Recombinant Mouse Monoclonal Antibody

Catalog # F104808

Product Information

Application	FC
Recommended Usage	5 μ L per million cells in 100 μ L staining volume or 5 μ L per 100 μ L of whole blood.
Reactivity	Human
Host	Mouse
Clonality	Monoclonal
Clone No.	11H73I47
Isotype	IgG2b, κ
Label	PE (Ex/Em: 565/574 nm)
Format	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and 0.2% (w/v) BSA.
Storage	Shipped on wet ice. Store undiluted between 2°C and 8°C and protected from prolonged exposure to light. Do not freeze.
Precautions	PE Mouse IgG2b, κ Isotype Ctrl Antibody [11H73I47] is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information
