

## PE Mouse IgG1, $\kappa$ Isotype Ctrl Antibody

Purified PE-conjugated Recombinant Mouse Monoclonal Antibody

Catalog # F104608

### Product Information

---

Application	FC
Recommended Usage	5 $\mu$ L per million cells in 100 $\mu$ L staining volume or 5 $\mu$ L per 100 $\mu$ L of whole blood.
Reactivity	Human
Host	Mouse
Clonality	Monoclonal
Clone No.	75S23A27
Isotype	IgG1, $\kappa$
Label	PE (Ex/Em: 565/574 nm)
Format	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and 0.2% (w/v) BSA.
Storage	Shipped on wet ice. Store undiluted between 2°C and 8°C and protected from prolonged exposure to light. Do not freeze.
Precautions	PE Mouse IgG1, $\kappa$ Isotype Ctrl Antibody [75S23A27] is for research use only and not for use in diagnostic or therapeutic procedures.

### Protein Information

---