

## FITC Mouse IgG2b, $\kappa$ Isotype Ctrl Antibody

Purified FITC-conjugated Recombinant Mouse Monoclonal Antibody

Catalog # F104805

### Product Information

---

Application	FC
Recommended Usage	5 $\mu$ L per million cells in 100 $\mu$ L staining volume or 5 $\mu$ L per 100 $\mu$ L of whole blood.
Reactivity	Human
Host	Mouse
Clonality	Monoclonal
Clone No.	11H73I47
Isotype	IgG2b, $\kappa$
Label	FITC (Ex/Em: 491/516 nm)
Format	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and 0.2% (w/v) BSA.
Storage	Shipped on wet ice. Store undiluted between 2°C and 8°C and protected from prolonged exposure to light. Do not freeze.
Precautions	FITC Mouse IgG2b, $\kappa$ Isotype Ctrl Antibody [11H73I47] is for research use only and not for use in diagnostic or therapeutic procedures.

### Protein Information

---