

VF450 Mouse IgG1, κ Isotype Ctrl Antibody

Purified VF450-conjugated Recombinant Mouse Monoclonal Antibody

Catalog # F104601

Product Information

Application	FC
Recommended Usage	5 μ L per million cells in 100 μ L staining volume or 5 μ L per 100 μ L of whole blood.
Reactivity	Human
Host	Mouse
Clonality	Monoclonal
Clone No.	75S23A27
Isotype	IgG1, κ
Label	VF450 (Ex/Em: 406/445 nm)
Format	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and 0.2% (w/v) BSA.
Storage	Shipped on wet ice. Store undiluted between 2°C and 8°C and protected from prolonged exposure to light. Do not freeze.
Precautions	VF450 Mouse IgG1, κ Isotype Ctrl Antibody [75S23A27] is for research use only and not for use in diagnostic or therapeutic procedures.

Protein Information