

## APC Mouse IgG2b, $\kappa$ Isotype Ctrl Antibody

Purified APC-conjugated Recombinant Mouse Monoclonal Antibody

Catalog # F104815

### Product Information

Application	FC
Recommended Usage	5 $\mu$ L per million cells in 100 $\mu$ L staining volume or 5 $\mu$ L per 100 $\mu$ L of whole blood.
Reactivity	Human
Host	Mouse
Clonality	Monoclonal
Clone No.	11H73147
Isotype	IgG2b, $\kappa$
Label	APC (Ex/Em: 651/660 nm)
Format	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and 0.2% (w/v) BSA.
Storage	Shipped on wet ice. Store undiluted between 2°C and 8°C and protected from prolonged exposure to light. Do not freeze.
Precautions	APC Mouse IgG2b, $\kappa$ Isotype Ctrl Antibody [11H73147] is for research use only and not for use in diagnostic or therapeutic procedures.

### Protein Information