

Anti-Nuclei Mouse mAb


Purified Recombinant Mouse Monoclonal Antibody

Catalog # M015800


Product Information

Application	ELISA
Reactivity	Human
Dilution	IHC-P 1:100~1:200; IF 1:100~1:200
Host	Mouse
Clonality	Monoclonal
Clone No.	84E15P56
Isotype	IgG
Label	Unconjugated
Immunogen	Purified nuclei from human cells.
Format	Affinity purified monoclonal antibody supplied in PBS with 0.02% sodium azide and 50% glycerol, pH 7.3.
Storage	Shipped on wet ice. Store at -20°C. Stable for 24 months from date of receipt. Aliquoting is unnecessary for -20°C storage.
Precautions	Anti-Nuclei Mouse mAb [84E15P56] is for research use only and not for use in diagnostic or therapeutic procedures.


Protein Information



Immunohistochemistry - Anti-Nuclei Mouse mAb [84E15P56]
 Sample: Paraformaldehyde-fixed, paraffin embedded human cervical cancer tissue
 Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0) for 30 mins.
 Primary antibody: M015800 at 1:200 dilution
 Secondary antibody: Goat Anti-Mouse IgG (H+L), HRP conjugated at 1:1,000 dilution
 DAB was used as the chromogen.
 Counter stained with hematoxylin.
 Positive/negative staining were presented.
 Only the secondary antibody was used as the negative control.



Immunohistochemistry - Anti-Nuclei Mouse mAb [84E15P56]
 Sample: Paraformaldehyde-fixed, paraffin embedded human endometrial carcinoma tissue
 Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0) for 30 mins.
 Primary antibody: M015800 at 1:200 dilution
 Secondary antibody: Goat Anti-Mouse IgG (H+L), HRP conjugated at 1:1,000 dilution
 DAB was used as the chromogen.
 Counter stained with hematoxylin.
 Positive/negative staining were presented.
 Only the secondary antibody was used as the negative control.



Immunofluorescence - Anti-Nuclei Mouse mAb [84E15P56]
 Sample: HeLa cells
 The cells were fixed with 4% paraformaldehyde (10 min), permeabilized with 0.5% Triton X-100 for 10 minutes and then blocked with 5% BSA in 0.1% PBS-Tween for 0.5 hours.
 Primary antibodies: M015800 at 1:100 dilution
 Secondary antibodies: Goat anti-Mouse (488) at 1:1,000 dilution (shown in green)
 Nuclei were stained with DAPI (shown in blue).