

| 货号 | 英文名称 | 货号 | 英文名称 |
|-------|-------------------------------------|-------|-------------------------------------|
| 小鼠 | | | |
| HJ153 | Mouse ACE ELISA Kit | HJ161 | Mouse EGF ELISA Kit |
| HJ176 | Mouse IL-1 α ELISA Kit | HJ177 | Mouse IL-1 β ELISA Kit |
| HJ178 | Mouse IL-2 ELISA Kit | HJ179 | Mouse IL-4 ELISA Kit |
| HJ181 | Mouse IL-5 ELISA Kit | HJ182 | Mouse IL-6 ELISA Kit |
| HJ185 | Mouse IL-10 ELISA Kit | HJ187 | Mouse IL-17A ELISA Kit |
| HJ196 | Mouse MCP-1 ELISA Kit | HJ197 | Mouse MMP-9 ELISA Kit |
| HJ205 | Mouse TGF- β 1 ELISA Kit | HJ207 | Mouse TNF- α ELISA Kit |
| HJ210 | Mouse VEGF ELISA Kit | HJ239 | Mouse BDNF ELISA Kit |
| HJ242 | Mouse IgG1 ELISA Kit | HJ245 | Mouse MPO ELISA Kit |
| HJ247 | Mouse IL-12p40 ELISA Kit | HJ248 | Mouse G-CSF ELISA Kit |
| HJ249 | Mouse GM-CSF ELISA Kit | HJ251 | Mouse IgG2a ELISA Kit |
| HJ252 | Mouse IgG2b ELISA Kit | HJ253 | Mouse IgG3 ELISA Kit |
| HJ254 | Mouse IL-13 ELISA Kit | HJ258 | Mouse CCL11 Eotaxin ELISA Kit |
| HJ257 | Mouse Angiopoietins/Ang-1 ELISA Kit | HJ261 | Mouse ALK-1/ACVRL1 ELISA Kit |
| HJ265 | Mouse β -NGF ELISA Kit | HJ268 | Mouse Activin A ELISA Kit |
| HJ270 | Mouse BMP-2 ELISA Kit | HJ273 | Mouse BMP-4 ELISA Kit |
| HJ275 | Mouse CDH3/P-cadherin ELISA Kit | HJ278 | Mouse IL-3 ELISA Kit |
| HJ280 | Mouse IL-20 ELISA Kit | HJ283 | Mouse CDH1/E-cadherin ELISA Kit |
| HJ191 | Mouse IL-23 ELISA Kit | HJ190 | Mouse IL-22 ELISA Kit |
| HJ192 | Mouse IL-33 ELISA Kit | HJ209 | Mouse VCAM-1 ELISA Kit |
| HJ287 | Mouse Angiopoietins/Ang-2 ELISA Kit | HJ288 | Mouse IL-15 ELISA Kit |
| HJ289 | Mouse IL-11 ELISA Kit | HJ290 | Mouse TNFRSF9/4-1BB/CD137 ELISA Kit |



扫码查看ELISA小贴士



上海雅酶生物医药科技有限公司
Shanghai EpiZyme Biomedical Technology Co., Ltd
地址:上海市闵行区新骏环路760号3层
网址:www.epizyme.cn 电话:400 058 8030

明星产品



ELISA检测试剂盒



01.

控制板间板内
差异 < 7%

02.

回收率范围
80-120%

03.

线性度范围
80-120%

04.

现货速达,
代测液体样本

严格的质检标准

产品原料行业金标准

稳定性高

检测样本种类多

准确率高

精密度高, 批间差小

重复性强, 回收率高

专业客服一对一指导


ELISA 实验对比数据

正确ELISA标准曲线的OD值范围一般在0.02-2.0之间，且呈现上升趋势。


| | |
|--|--|
| <p>雅酶的ELISA标准曲线：</p> <p>符合正确ELISA标准曲线范围且呈现上升的趋势</p> | <p>不理想的ELISA标准曲线：</p> <p>上升趋势不明显，且该标准曲线误差较大</p> |
|--|--|

Mouse IL-6 ELISA KIT

雅酶的ELISA标准曲线数据：




不理想的ELISA标准曲线数据：




▲ 上升趋势不明显，且最高OD值没有到达2左右

Mouse IL-10 ELISA KIT

雅酶的ELISA标准曲线数据：




不理想的ELISA标准曲线数据：




▲ 上升趋势不明显，且标准曲线误差较大

Human VEGF ELISA KIT

雅酶的ELISA标准曲线数据：



不理想的ELISA标准曲线数据：



▲ 上升趋势不明显，且标准曲线误差较大

ELISA 部分精品展示



查看更多
请扫二维码

| 货号 | 英文名称 | 货号 | 英文名称 |
|-------|-------------------------------------|-------|-------------------------------------|
| HJ004 | Human Angiopoietins/Ang-1 ELISA Kit | HJ005 | Human Angiopoietins/Ang-2 ELISA Kit |
| HJ009 | Human BMP-2 ELISA Kit | HJ018 | Human CCL11/Eotaxin ELISA Kit |
| HJ031 | Human EGF ELISA Kit | HJ039 | Human G-CSF ELISA Kit |
| HJ042 | Human GM-CSF ELISA Kit | HJ048 | Human IFN-γ ELISA Kit |
| HJ059 | Human IL-1α ELISA Kit | HJ060 | Human IL-1β ELISA Kit |
| HJ061 | Human IL-2 ELISA Kit | HJ062 | Human IL-4 ELISA Kit |
| HJ063 | Human IL-5 ELISA Kit | HJ064 | Human IL-6 ELISA Kit |
| HJ066 | Human IL-8 ELISA Kit | HJ071 | Human IL-11 ELISA Kit |
| HJ072 | Human IL-12p70 ELISA Kit | HJ074 | Human IL-17 ELISA Kit |
| HJ078 | Human IL-23 ELISA Kit | HJ079 | Human IL-33 ELISA Kit |
| HJ085 | Human MCP-1 ELISA KIT | HJ090 | Human MMP-3 ELISA Kit |
| HJ091 | Human MMP-9 ELISA Kit | HJ092 | Human MPO ELISA Kit |
| HJ094 | Human β-NGF ELISA Kit | HJ097 | Human PCT ELISA Kit |

| | | | |
|-------|-------------------------------------|-------|-------------------------------------|
| HJ106 | Human TGF-β1 ELISA Kit | HJ110 | Human TNF-α ELISA Kit |
| HJ114 | Human VCAM-1 ELISA Kit | HJ115 | Human VEGF ELISA Kit |
| HJ119 | Human Insulin ELISA Kit | HJ120 | Human ACE ELISA Kit |
| HJ230 | Human NT-proBNP ELISA Kit | HJ231 | Human NGAL ELISA Kit |
| HJ232 | Human Lp-PLA2 ELISA Kit | HJ233 | Human L-FABP ELISA Kit |
| HJ234 | Human ProGRP ELISA Kit | HJ235 | Human KIM-1/TIM-1/HAVCR-1 ELISA Kit |
| HJ236 | Human GDF15 ELISA Kit | HJ237 | Human BDNF ELISA Kit |
| HJ240 | Human IL-18 ELISA Kit | HJ241 | Human IL-7 ELISA Kit |
| HJ243 | Human FGF-2 ELISA Kit | HJ244 | Human Serpin G1 ELISA Kit |
| HJ246 | Human COL1A1 ELISA Kit | HJ260 | Human ALK-1/ACVRL1 ELISA Kit |
| HJ262 | Human JAML/AMICA ELISA Kit | HJ263 | Human Amphiregulin/AREG ELISA Kit |
| HJ264 | Human ANGPTL4 ELISA Kit | HJ266 | Human Activin A ELISA Kit |
| HJ271 | Human BMP-4 ELISA Kit | HJ274 | Human CDH3/P-cadherin ELISA Kit |
| HJ276 | Human IL-3 ELISA Kit | HJ279 | Human IL-20 ELISA Kit |
| HJ281 | Human TNFRSF9/4-1BB/CD137 ELISA Kit | HJ282 | Human IL-16 ELISA Kit |
| HJ286 | Human IL-15 ELISA Kit | HJ291 | Human NKR2B4/CD244 ELISA Kit |
| 大鼠 | | | |
| HJ130 | Rat IFN-γ ELISA Kit | HJ134 | Rat IL-1α ELISA Kit |
| HJ135 | Rat IL-1β ELISA Kit | HJ136 | Rat IL-2 ELISA Kit |
| HJ139 | Rat IL-10 ELISA Kit | HJ141 | Rat IL-18 ELISA Kit |
| HJ144 | Rat MCP-1 ELISA Ki | HJ147 | Rat TGF-β1 ELISA Kit |
| HJ149 | Rat TNF-α ELISA Kit | HJ150 | Rat VEGF ELISA Kit |
| HJ238 | Rat BDNF ELISA Kit | HJ255 | Rat EGF ELISA Kit |
| HJ256 | Rat IL-13 ELISA Kit | HJ259 | Rat β-NGF ELISA Kit |
| HJ267 | Rat Activin A ELISA Kit | HJ269 | Rat BMP-2 ELISA Kit |
| HJ272 | Rat BMP-4 ELISA Kit | HJ277 | Rat IL-3 ELISA Kit |
| HJ284 | Rat GM-CSF ELISA Kit | HJ285 | Rat IL-5 ELISA Kit |