

雅酶精品抗体

实事求是做好每一次验证实验

重组兔单克隆抗体平台

噬菌体展示平台

鼠杂交瘤单抗平台

兔多克隆抗体平台

科学严谨的检测体系 专业的技术支持团队



高特异性






高亲和力



高批间稳定性

雅酶精品抗体结果展示

Anti-Phospho-Histone H2A.X (Ser139) Rabbit mAb [15L14K33]



雅酶公众号



雅酶视频号

📍 地址：上海市闵行区新骏环路760号1号楼3楼
🌐 网址：www.epizyme.cn
☎ 电话：400-058-8030
✉ 邮箱：info@epizyme.cn

部分热卖抗体

货号	名称	应用	种属
R013789	Anti-Phospho-NF-kB p65 (Ser536) Rabbit mAb [12L74M11]	WB,IHC-P/IF (Tissue-P),IF (Cell)/ICC,ELISA	Mouse,Rat,Human
R012746	Anti-NF-KB p65 Rabbit mAb [99M06L21]	WB,IF (Cell)/ICC,IHC-P/IF (Tissue-P),ELISA	Mouse,Human
R011234	Anti-AKT1/2/3 Rabbit mAb [69K96M86]	WB,IHC-P/IF (Tissue-P),ELISA	Human,Mouse
R013644	Anti-Bax Rabbit mAb [59M52K38]	WB,IHC-P/IF (Tissue-P),ELISA	Human,Mouse,Rat
R011716	Anti-Glutathione Peroxidase 4 Rabbit mAb [82K58M44]	WB,IHC-P/IF (Tissue-P),ELISA	Human,Mouse,Rat
R013822	Anti-alpha smooth muscle Actin Rabbit mAb [30L29F93]	ELISA,WB,IHC-P/IF (Tissue-P),IF (Cell)/ICC	Human,Mouse,Rat
R013940	Anti-LC3B Rabbit mAb [95F45F38]	ELISA,WB,IHC-P/IF (Tissue-P),IF (Cell)/ICC	Mouse,Rat,Human
R013501	Anti-ERK1/2 Rabbit mAb [13M15K92]	WB,IHC-P/IF (Tissue-P),ELISA	Rat,Mouse,Human
R013735	Anti-PD-L1 Rabbit mAb [31G76K98]	ELISA,WB,IHC-P/IF (Tissue-P)	Human,Mouse
R011470	Anti-Phospho-AKT (Ser473) Rabbit mAb [81M88K84]	IHC-P/IF (Tissue-P),ELISA,WB	Human,Mouse
R013790	Anti-SQSTM1/p62 Rabbit mAb [53N86E71]	ELISA,WB,IHC-P/IF (Tissue-P)	Rat,Human,Mouse
R010101	Anti-Phospho-Histone H2A.X (Ser139) Rabbit mAb [15L14K33]	ELISA,IF (Cell)/ICC,WB,IHC-P/IF (Tissue-P)	Mouse,Rat,Human
R900001	Anti-BCL2 Rabbit mAb [47A21F24]	WB,ELISA	Human
R013743	Anti-Ki67 Rabbit mAb [36E24F01]	IHC-P/IF (Tissue-P),IF (Cell)/ICC,ELISA	Human,Mouse,Rat
R010732	Anti-Phospho-mTOR (Ser2448) Rabbit mAb [81K48K13]	IHC-P/IF (Tissue-P),IF (Cell)/ICC,ELISA,WB	Rat,Human,Mouse
R013244	Anti-E Cadherin Rabbit mAb [10M94M96]	WB,IHC-P/IF (Tissue-P),ELISA,IF (Cell)/ICC	Human
R013028	Anti-PI3 Kinase p85 alpha Rabbit mAb [78L02M83]	WB,IHC-P/IF (Tissue-P),ELISA	Mouse,Rat,Human
R014250	Anti-Lamin B1 Rabbit mAb [54I84Q27]	WB,IHC-P/IF (Tissue-P),IF (Cell)/ICC,ELISA	Human,Mouse,Rat
R013762	Anti-Smad2/3 Rabbit mAb [14L86C22]	ELISA,IF (Cell)/ICC,WB,IHC-P/IF (Tissue-P)	Human,Mouse,Rat
R015076	Anti-HIF-1 alpha Rabbit mAb [39H77C67]	WB,IHC-P/IF (Tissue-P),IF (Cell)/ICC,ELISA	Human,Mouse,Rat
R011421	Anti-N Cadherin Rabbit mAb [14L59L66]	ELISA,IHC-P/IF (Tissue-P),WB	Human,Mouse,Rat
M010697	Anti-Vimentin Mouse mAb [19M28K95]	ELISA,IF (Cell)/ICC,WB,IHC-P/IF (Tissue-P)	Human
R013285	Anti-FACL4 Rabbit mAb [99M02K23]	IF (Cell)/ICC,ELISA,WB,IHC-P/IF (Tissue-P)	Rat,Human,Mouse
R010292	Anti-beta Catenin Rabbit mAb [41L48K19]	WB,IF (Cell)/ICC,ELISA,IHC-P/IF (Tissue-P)	Human,Mouse,Rat
R010272	Anti-JNK1/2/3 Rabbit mAb [62M83L63]	WB,IHC-P/IF (Tissue-P),ELISA	Human,Mouse,Rat
R013794	Anti-p38 Rabbit mAb [49C01B92]	ELISA,WB,IHC-P/IF (Tissue-P),IF (Cell)/ICC	Human,Mouse,Rat
R012092	Anti-Caspase 3 Rabbit mAb [32M89K63]	WB,ELISA,IHC-P/IF (Tissue-P)	Human
R011464	Anti-TSG101 Rabbit mAb [22L57K84]	IF (Cell)/ICC,ELISA,WB	Human,Mouse,Rat
R013921	Anti-CD9 Rabbit mAb [85E05T23]	WB,ELISA	Rat,Mouse,Human
R900002	Anti-ZO-1 Rabbit mAb [88E68I39]	ELISA,WB,IHC-P/IF (Tissue-P),IF (Cell)/ICC	Human,Mouse,Rat
R014764	Anti-Heme Oxygenase 1 Rabbit mAb [58Q07H13]	ELISA,WB,IHC-P/IF (Tissue-P)	Mouse,Rat,Human